



SHOULDER STABILIZATION SYSTEM Measuring Guide



MEASURING:

1. Measure from the Point of Elbow on one side, over the back, to Point of Elbow on the opposite side.
2. Measure around the leg at the Point of Elbow.
3. Measure around the leg 4 inches below the Point of Elbow.
4. Measure from the Point of Elbow to the Top of Paw on either leg.

** Alternate Measurement #3 (Both Legs) - If Measurement #4 is 6 inches or less, measure the leg half way between the Point of Elbow and Top of Paw

